

Product datasheet for **RC203085A1V**

Human DUSP13 (NM_016364) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human DUSP13 (NM_016364) AAV Particle
Tag: Myc-DDK
Symbol: DUSP13
Synonyms: BEDP; DUSP13A; DUSP13B; MDSP; SKRP4; TMDP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203085 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGACTCACTGCAGAAGCAGGACCTCCGGAGGCCAAGATCCATGGGGCAGTCCAGGCATCTCCCTACC
AGCCGCCACATTGGCTTCGCTGCAGCGCTTGTGTGGGTCCGTCAAGCTGCCACACTGAACCATATCGA
TGAGGTCTGGCCAGCCTCTTCTGGGAGATGCGTACGCAGCCCGGACAAGAGCAAGCTGATCCAGCTG
GGAATCACCCACGTTGTGAATGCCGCTGCAGGCAAGTTCAGGTGGACACAGGTGCCAAATCTACCGTG
GAATGTCCCTGGAGTACTATGGCATCGAGGCGGATGACAACCCCTTCTTCGACCTCAGTGTCTACTTTCT
GCCTGTTGCTCGATACATCCGAGCTGCCCTCAGTGTTCCTCAAGGCCGCGTGTGTACACTGTGCCATG
GGGGTAAGCCGCTCTGCCACACTTGTCTGGCCTTCTCATGATCTATGAGAACATGACGCTGGTAGAGG
CCATCCAGACGGTGCAGGCCACCGCAATATCTGCCCTAACTCAGGCTTCTCCGGCAGCTCCAGGTTCT
GGACAACCGACTGGGGCGGGAGACGGGGCGGTTTC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203085 protein sequence
 Red=Cloning site Green=Tags(s)

MDSLQKQDLRRPKIHGAVQASPYQPPTLASLQRLWVRQAATLNHIDEVWPSLFLGDAYAARDKSKLIQL
GITHVNAAGKFQVDTGAKFYRGMSLEYGIEADNPFFDLVYFLPVARYIRAALSVPGRVLVHCAM
GVSRATLVLAFLMIYENMTLVEAIQTVQAHNRNICPNSGFLRQLQVLDNRLGRETGRF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human


[View online »](#)

Serotype:	AAV-2
ACCN:	NM_016364
ORF Size:	594 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_016364.3</u>
RefSeq Size:	923 bp
RefSeq ORF:	597 bp
Locus ID:	51207
UniProt ID:	<u>Q9UII6</u>
Cytogenetics:	10q22.2
MW:	22.2 kDa