

Product datasheet for **RC203099A1V**

Human ZNF193 (ZSCAN9) (NM_006299) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ZNF193 (ZSCAN9) (NM_006299) AAV Particle
Tag:	Myc-DDK
Symbol:	ZNF193
Synonyms:	PRD51; ZNF193
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203099 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAATACAACTCAAAGGAGGTTTTATCCCTGGGTGTTCAAGTCCCGAGGCATGGGAAGAACTTCTGA
CAATGAAAGTGAAGCAAAAAGTCACCTTCAATGGCAGGAATCCAGACTGAAACGCAGTAATCCACTGGC
AAGGGAATCTTCCGAAGGCACCTTCGACAGCTGTGCTACCAAGAGACCCCTGGACCAAGGGAGGCTCTT
ACTCGACTCCAGGAACCTTGTACCAAGTGGTTGAGGCCACATGTGAGCACAAAGGAGCAGATTTTGGATC
TGCTGGTGTGGAGCAGTTTCTATCCATTCTGCCAAGGAGCTCCAGGGCTGGGTGAGGGAACACTGTCC
AGAGAGTGGAGAAGAGGCTGTGATTTTCTGGAGGATCTGGAGAGAGAGCTCGATGAACCACAACATGAG
ATGGTGGCCACAGACACAGACAAGAAGTCCTCTGTAAAGAGATGGTGCCTCTAGCAGAGCAGACACCAC
TGACCCTTCAGTCCCAGCCTAAGGAGCCACAGCTCACATGTGACTCTGCTCAGAAGTGCCATTCTATTGG
AGAGACAGATGAAGTAACCAAGACTGAGGACAGAGAGTTGGTGCTAAGGAAAGACTGTCCTAAGATAGTG
GAACCACATGGGAAAATGTTAATGAGCAGACCTGGGAGGTATCACAGCAGGATCCCTCACATGGAGAAG
TTGGTGAACATAAGGATAGGATAGAGAGGCAGTGGGAAACCTCTTAGGAGAGGGGAACACAAATGTGA
TGAATGTGGGAAGAGCTTACTCAGAGCTCAGGTCTCATTGACATCAAAGAATTCATACTGGAGAAAGA
CCTTATGAATGTAATGAATGTGGAAAGCCTTCAGTCGAAGTTCTGGTCTTTTAAATCACCGAGGAATCC
ACAATATACAGAAACGGTACCACTGCAAGGAGTGTGGGAAGGTCTTCAGTCAGAGTGCGGGTCTTATCCA
GCATCAGAGAATCCACAAGGAGAAAAGCCGTATCAGTCAGCCAGTGCAGTAAGAGCTACAGTCGGCGT
TCATTTCTCATTGAACATCAGAGAAGCCACAGGGGAGCGACCTCACCAGTGCATTGAATGTGGGAAAA
GCTTTAATCGACACTGCAACCTCATTGCCATCAGAAGATCCACACAGTGGCTGAGCTGGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203099 protein sequence Red=Cloning site Green=Tags(s) MNTNSKEVLSLGVQVPEAWEELLTMKVEAKSHLQWQESRLKRSNPLAREIFRRHFRQLCYQETPGPREAL TRLQELCYQWLRPHVSTKEQILDLLVLEQFLSILPKELQGWVREHCPESGEEAVILLEDLERELDEPQHE MVAHRHRQEVLCHEMVPLAEQTPLTLQSQPKPQLTCDSAQKCHSIGETDEVTKTEDRELVLKDCPKIV EPHGKMFNEQTWEVSQQDP SHGEVGEHKDRIERQWGNLLGEGQHKCECGKSFTQSSGLIRHQRIHTGER PYECNECGKAFSRSSGLFNHRGIHNIQKRYHCKEKGKVFQSAGLIQHQRHKGKPYQCSQCSKSYSRR SFLIEHQRSHTGERPHQCIIECGKSFNRHCNLIHQKIHTVAELV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006299
ORF Size:	1182 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006299.2
RefSeq Size:	1698 bp
RefSeq ORF:	1185 bp
Locus ID:	7746
UniProt ID:	O15535
Cytogenetics:	6p22.1
MW:	46 kDa