

## Product datasheet for **RC203117L2V**

### Nuclear Matrix Protein p84 (THOC1) (NM\_005131) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Nuclear Matrix Protein p84 (THOC1) (NM_005131) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Nuclear Matrix Protein p84
Synonyms:	HPR1; P84; P84N5
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_005131
ORF Size:	1971 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203117).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005131.1</a>
RefSeq Size:	2133 bp
RefSeq ORF:	1974 bp
Locus ID:	9984
UniProt ID:	<a href="#">Q96FV9</a>
Cytogenetics:	18p11.32
Domains:	DEATH



[View online »](#)

**Protein Pathways:** Proteasome, Spliceosome

**MW:** 75.7 kDa

**Gene Summary:** HPR1 is part of the TREX (transcription/export) complex, which includes TEX1 (MIM 606929), THO2 (MIM 300395), ALY (MIM 604171), and UAP56 (MIM 142560).[supplied by OMIM, Nov 2010]