

Product datasheet for **RC203123L4V**

HTN3 (NM_000200) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HTN3 (NM_000200) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HTN3
Synonyms:	HIS2; HTN2; HTN5; PB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000200
ORF Size:	153 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203123).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000200.1
RefSeq Size:	601 bp
RefSeq ORF:	156 bp
Locus ID:	3347
UniProt ID:	P15516
Cytogenetics:	4q13.3
Protein Families:	Secreted Protein
MW:	6.1 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the histatin family of small, histidine-rich, cationic proteins. They function as antimicrobial peptides and are important components of the innate immune system. Histatins are found in saliva and exhibit antibacterial, antifungal activities and function in wound healing. [provided by RefSeq, Sep 2014]