

Product datasheet for **RC203125L3V**

KCNJ15 (NM_002243) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KCNJ15 (NM_002243) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KCNJ15
Synonyms:	IRKK; KIR1.3; KIR4.2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002243
ORF Size:	1125 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203125).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002243.3
RefSeq Size:	2942 bp
RefSeq ORF:	1128 bp
Locus ID:	3772
UniProt ID:	Q99712
Cytogenetics:	21q22.13-q22.2
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane
MW:	42.6 kDa



[View online »](#)

Gene Summary:

Potassium channels are present in most mammalian cells, where they participate in a wide range of physiologic responses. The protein encoded by this gene is an integral membrane protein and inward-rectifier type potassium channel. The encoded protein has a greater tendency to allow potassium to flow into a cell rather than out of a cell. Eight transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Feb 2013]