

Product datasheet for **RC203128A1V**

Human Tetranectin (CLEC3B) (NM_003278) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Tetranectin (CLEC3B) (NM_003278) AAV Particle
Tag: Myc-DDK
Symbol: Tetranectin
Synonyms: TN; TNA
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203128 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGCTCTGGGGGGCCTACCTCCTCCTGCCTCTTCTCCCTCCTGACCCAGGTACCAACCGAGCCAC
 CAACCCAGAAGCCCAAGAAGATTGTAATGCCAAGAAAGATGTTGTGAACACAAAGATGTTTGAGGAGCT
 CAAGAGCCGCTGGACACCCTGGCCAGGAGGTGGCCCTGCTGAAGGAGCAGCAGGCCCTGCAGACGGTC
 TGCCTGAAGGGGACCAAGGTGCATGAAATGCTTTCTGGCCTTCACCCAGACGAAGACCTTCACGAGG
 CCAGCGAGGACTGCATCTCGCGCGGGGGCACCCTGAGCACCCCTCAGACTGGCTCGGAGAACGACGCCCT
 GTATGAGTACCTGCGCCAGAGCGTGGGCAACGAGGCCGAGATCTGGCTGGGCCTCAACGACATGGCGGCC
 GAGGGCACCTGGGTGGACATGACCGGCGCCCGCATCGCCTACAAGAACTGGGAGACTGAGATCACCGCGC
 AACCCGATGGCGGCAAGACCGAGAACTGCGCGGTCTGTGAGGCGCGGCAACGGCAAGTGGTTCGACAA
 GCGCTGCCGCGATCAGCTGCCCTACATCTGCCAGTTCGGGATCGTG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203128 protein sequence
 Red=Cloning site Green=Tags(s)

MELWGAYLLLCLFSLLTQVTTEPPTQPKKIVNAKDVVNTKMFEELKSRLDTLAQEVALLEQQALQTV
 CLKGTQVHMKFLAFTQTKTFHEASEDCISRGTLSTPQTGSENDALYEYLRQSVGNEAEIWLGLNDMAA
 EGTWVDMTGARIAYKNWETEITAQPDGGKTENCAVLGAANGKWFDRKCRDQLPYICQFGIV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_003278
ORF Size:	606 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003278.1</u> , <u>NP_003269.1</u>
RefSeq Size:	886 bp
RefSeq ORF:	609 bp
Locus ID:	7123
UniProt ID:	<u>P05452</u>
Cytogenetics:	3p21.31
MW:	22.6 kDa