

Product datasheet for **RC203130L3V**

SR1 (SRI) (NM_003130) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SR1 (SRI) (NM_003130) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SR1 |
| Synonyms: | CP-22; CP22; SCN; V19 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_003130 |
| ORF Size: | 594 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC203130). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_003130.2 |
| RefSeq Size: | 2036 bp |
| RefSeq ORF: | 597 bp |
| Locus ID: | 6717 |
| UniProt ID: | P30626 |
| Cytogenetics: | 7q21.12 |
| Domains: | EFh |
| Protein Families: | Druggable Genome |



[View online »](#)

MW: 21.7 kDa

Gene Summary: This gene encodes a calcium-binding protein with multiple E-F hand domains that relocates from the cytoplasm to the sarcoplasmic reticulum in response to elevated calcium levels. In addition to regulating intracellular calcium homeostasis it also modulates excitation-contraction coupling in the heart. Alternative splicing results in multiple transcript variants encoding distinct proteins. Multiple pseudogenes exist for this gene. [provided by RefSeq, Mar 2012]