

## Product datasheet for **RC203135L3V**

### GCDFP 15 (PIP) (NM\_002652) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GCDFP 15 (PIP) (NM_002652) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GCDFP 15
Synonyms:	GCDFP-15; GCDFP15; GPIP4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002652
ORF Size:	438 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203135).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002652.2</a>
RefSeq Size:	591 bp
RefSeq ORF:	441 bp
Locus ID:	5304
UniProt ID:	<a href="#">P12273</a>
Cytogenetics:	7q34
Protein Families:	Secreted Protein
MW:	16.6 kDa



[View online »](#)