

Product datasheet for **RC203142A1V**

Human Tissue Factor (F3) (NM_001993) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Tissue Factor (F3) (NM_001993) AAV Particle
Tag:	Myc-DDK
Symbol:	Tissue Factor
Synonyms:	CD142; TF; TFA
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203142 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGACCCCTGCCTGGCCCCGGGTCCCGCGCCCCGAGACCGCCGTCGCTCGGACGCTCCTGCTCGGCT
GGGTCTTCGCCCAGGTGGCCGGCGCTTCAGGCACTACAAATACTGTGGCAGCATATAATTTAACTTGAA
ATCAACTAATTTCAAGACAATTTGGAGTGGGAACCAACCCGTCAATCAAGTCTACACTGTTCAAATA
AGCACTAAGTCAGGAGATTGAAAAGCAAATGCTTTTACACAACAGACACAGAGTGTGACCTACCGGACG
AGATTGTGAAGGATGTGAAGCAGACGTACTTGGCAGGGTCTTCTCCTACCCGGCAGGGAATGTGGAGAG
CACCGGTTCTGCTGGGGAGCCTCTGTATGAGAACTCCCCAGAGTTCACACCTTACCTGGAGACAAACCTC
GGACAGCCAACAATTCAGAGTTTTGAACAGGTGGGAACAAAAGTGAATGTGACCGTAGAAGATGAACGGA
CTTTAGTCAGAAGGAACAACACTTTCTAAGCCTCCGGGATGTTTTTGGCAAGGACTTAATTTATACACT
TTATTATTGGAATCTTCAAGTTCAGGAAAGAAAACAGCCAAAACAACTAATGAGTTTTTGATTGAT
GTGGATAAAGGAGAAACTACTGTTTCAGTGTTCAGCAGTGATTCCCTCCCGAACAGTTAACCGGAAGA
GTACAGACAGCCCGTAGAGTGTATGGCCAGGAGAAAGGGAATTCAGAGAAATATTCTACATCATTGG
AGCTGTGGTATTTGTGGTCATCATCCTTGTATCATCCTGGCTATATCTCTACACAAGTGTAGAAAGGCA
GGAGTGGGGCAGAGCTGGAAGGAGAACTCCCCACTGAATGTTTCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203142 protein sequence Red=Cloning site Green=Tags(s) METPAWPRVPRPETAVARTLLLGWVFAQVAGASGTTNTVAAYNLTKSTNFKTILEWEPKPVNQVYTVQI STKSGDWKSKCFYTTDTECDLTDEIVKDVKQTYLARVFSYPAGNVESTGSAGEPLYENSPEFTPYLETNL GQPTIQSFEQVGTKVNVTVEDERTLVRRNNTFLSLRDVFGKDLIYTLYYWKSSSSGKKTAKTNTNEFLID VDKGENYCFSVQAVIPSRTVNRKSTDSPVECMGQEKGEFREIFYIIGAVVFVVIILVIIAISLHKCRKA GVGQSWKENSPLNVS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001993
ORF Size:	885 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001993.2
RefSeq Size:	2393 bp
RefSeq ORF:	888 bp
Locus ID:	2152
UniProt ID:	P13726
Cytogenetics:	1p21.3
MW:	33.1 kDa