

#### OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for RC203142L4

## Tissue Factor (F3) (NM\_001993) Human Tagged Lenti ORF Clone

## **Product data:**

Product Type:	Expression Plasmids
Product Name:	Tissue Factor (F3) (NM_001993) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Tissue Factor
Synonyms:	CD142; TF; TFA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203142).
<b>Restriction Sites:</b>	Sgfl-Mlul
Cloning Scheme:	
	Cloning sites used for ORF Shuttling:
	Sgf I         ORF         Mlu I            GCG ATC GC         ATG//         NNN         ACG CGT



ACCN: ORF Size: NM\_001993 885 bp



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


Service Pactor (F3) (NM_001993) Human Tagged Lenti ORF Clone – RC203142L4	
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol> <li>Centrifuge at 5,000xg for 5min.</li> <li>Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>Close the tube and incubate for 10 minutes at room temperature.</li> <li>Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
RefSeq:	<u>NM 001993.2</u>
RefSeq Size:	2393 bp
RefSeq ORF:	888 bp
Locus ID:	2152
UniProt ID:	<u>P13726</u>
Cytogenetics:	1p21.3
Domains:	Tissue_fac
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Complement and coagulation cascades
MW:	33.1 kDa

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

#### Scrigene Tissue Factor (F3) (NM\_001993) Human Tagged Lenti ORF Clone – RC203142L4

Gene Summary:This gene encodes coagulation factor III which is a cell surface glycoprotein. This factor<br/>enables cells to initiate the blood coagulation cascades, and it functions as the high-affinity<br/>receptor for the coagulation factor VII. The resulting complex provides a catalytic event that is<br/>responsible for initiation of the coagulation protease cascades by specific limited proteolysis.<br/>Unlike the other cofactors of these protease cascades, which circulate as nonfunctional<br/>precursors, this factor is a potent initiator that is fully functional when expressed on cell<br/>surfaces, for example, on monocytes. There are 3 distinct domains of this factor: extracellular,<br/>transmembrane, and cytoplasmic. Platelets and monocytes have been shown to express this<br/>coagulation factor III has been described. Platelet-dependent monocyte expression of<br/>coagulation factor III has been described to be associated with Coronavirus Disease 2019<br/>(COVID-19) severity and mortality. This protein is the only one in the coagulation pathway for<br/>which a congenital deficiency has not been described. Alternate splicing results in multiple<br/>transcript variants.[provided by RefSeq, Aug 2020]

## **Product images:**



Circular map for RC203142L4

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US