

Product datasheet for RC203146L4

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

hCG (CGA) (NM_000735) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: hCG (CGA) (NM_000735) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: hCG

Synonyms: CG-ALPHA; FSHA; GPA1; GPHa; GPHA1; HCG; LHA; TSHA

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC203146).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_000735

ORF Size: 348 bp



hCG (CGA) (NM_000735) Human Tagged Lenti ORF Clone - RC203146L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000735.2</u>

RefSeq Size: 768 bp
RefSeq ORF: 351 bp
Locus ID: 1081

UniProt ID: P01215
Cytogenetics: 6q14.3

Domains: GHA

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein

Protein Pathways: Autoimmune thyroid disease, GnRH signaling pathway, Neuroactive ligand-receptor

interaction

MW: 13 kDa

Gene Summary: The four human glycoprotein hormones chorionic gonadotropin (CG), luteinizing hormone


(LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) are dimers consisting of alpha and beta subunits that are associated noncovalently. The alpha subunits of these hormones are identical, however, their beta chains are unique and confer biological

specificity. The protein encoded by this gene is the alpha subunit and belongs to the glycoprotein hormones alpha chain family. Two transcript variants encoding different

isoforms have been found for this gene. [provided by RefSeq, Nov 2011]



Product images:



Circular map for RC203146L4