

## Product datasheet for **RC203148L3V**

### **BHMT (NM\_001713) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BHMT (NM_001713) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BHMT
Synonyms:	BHMT1; HEL-S-61p
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001713
ORF Size:	1218 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203148).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001713.1</a>
RefSeq Size:	2515 bp
RefSeq ORF:	1221 bp
Locus ID:	635
UniProt ID:	<a href="#">Q93088</a>
Cytogenetics:	5q14.1
Domains:	S-methyl_trans



[View online »](#)

<b>Protein Pathways:</b>	Cysteine and methionine metabolism, Glycine, serine and threonine metabolism, Metabolic pathways
<b>MW:</b>	45 kDa
<b>Gene Summary:</b>	This gene encodes a cytosolic enzyme that catalyzes the conversion of betaine and homocysteine to dimethylglycine and methionine, respectively. Defects in this gene could lead to hyperhomocyst(e)inemia, but such a defect has not yet been observed. [provided by RefSeq, Jul 2008]