

Product datasheet for **RC203152A1V**

Human trfp (MED20) (NM_004275) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human trfp (MED20) (NM_004275) AAV Particle
Tag: Myc-DDK
Symbol: trfp
Synonyms: PRO0213; SRB2; TRFP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203152 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGAGTGACTTGTGTGTCCAGATGCCTGTGGCCGAGGGCAAGAGTGTTAGCAAACCGTAGAGCTCC
 TTACCCGGAAATTGGAGATGCTTGGGGCAGAGAAGCAAGGAACATTTGTGTGGACTGTGAGACTTACCA
 TACGGCCGCTCTACCTTGGCAGCCAAGGTCAGACCGGGAAGCTGATGTATGTGATGCACAACCTCAGAG
 TACCCATTGAGCTGTTTCGCCCTCTTTGAGAATGGCCCTTGCCCTATTGCTGACACCAACTTTGATGTGC
 TTATGGTGAAGCTCAAGGGCTTTTCCAGAGTGCTAAGGCCAGCAAGATTGAGACCCGGGGCACCAGGTA
 CCAGTACTGTGACTTCTGTTGAAGGTGGGCACGGTCACAATGGGGCCAGTGCCCGGGGCATCTCTGTG
 GAGGTGGAGTATGGCCCTGTGTGGTAGCTAGTGACTGCTGGAGTCTGCTGCTCGAGTTCCTACAGAGTT
 TTCTAGGCAGCCACACACCAGGGGCTCCCGCAGTGTTTGGGAACAGACATGATGCGGTCTACGGCCAGC
 AGATACCATGGTCCAGTACATGGAACCTTTCAACAAGATCCGCAAGCAGCAGCAGGTGCCGGTGGCTGG
 ATTCGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203152 protein sequence
 Red=Cloning site Green=Tags(s)

MGVTCVSQMPVAEGKSVQQTVELLTKLEMLGAEKQGTFCVDCETYHTAASTLGSQGQTGKLMYVMHNS
 YPLSCFALFENGPCLIADTNFDVLMVKLGFFQSAKASKIETRTRYQYCDLVKVGTVTMGPSARGISV
 EVEYGPCVVASDCWSLLLEFLQSFLGSHTPGAPAVFGNRHDAVYGPADTMVQYMELFNKIRKQQQVPVAG
 IR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_004275
ORF Size:	636 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004275.3
RefSeq Size:	2478 bp
RefSeq ORF:	639 bp
Locus ID:	9477
UniProt ID:	Q9H944
Cytogenetics:	6p21.1
MW:	23.2 kDa