

## Product datasheet for **RC203170L2V**

### **SPIB (NM\_003121) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SPIB (NM_003121) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPIB
Synonyms:	SPI-B
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_003121
ORF Size:	786 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203170).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003121.2</a> , <a href="#">NP_003112.2</a>
RefSeq Size:	3537 bp
RefSeq ORF:	789 bp
Locus ID:	6689
UniProt ID:	<a href="#">Q01892</a>
Cytogenetics:	19q13.33
Protein Families:	Transcription Factors
MW:	28.8 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a transcriptional activator that binds to the PU-box (5'-GAGGAA-3') and acts as a lymphoid-specific enhancer. Four transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]