

Product datasheet for **RC203178A1V**

Human ISG20 (NM_002201) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Human ISG20 (NM_002201) AAV Particle
 Tag: Myc-DDK
 Symbol: ISG20
 Synonyms: CD25; HEM45
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >RC203178 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTGGGAGCCGTGAGGTGGTGGCCATGGACTGCGAGATGGTGGGGCTGGGGCCCCACCGGGAGAGTG
 GCCTGGCTCGTTGCAGCCTCGTGAACGTCCACGGTGTGTGTGTACGACAAGTTATCCGGCTGAGGG
 AGAGATACCGATTACAGAACCCGGGTGAGCGGGTCAACCCCTCAGCACATGGTGGGGGCCACACCATTT
 GCCGTGGCCAGGCTAGAGATCCTGCAGTCTCTGAAAGGCAAGCTGGTGGTGGGTGATGACCTGAAGCAGC
 ACTTCCAGGCACTGAAAGAGGACATGAGCGGCTACACAATCTACGACACGTCCACTGACAGGCTGTTGTG
 GCGTGAGGCCAAGCTGGACCACTGCAGGCGTGTCTCCCTGCGGGTGTGAGTGAGCGCCTCTGCACAAG
 AGCATCCAGAACAGCCTGCTTGGACACAGCTCGGTGGAAGATGCGAGGGCAACGATGGAGCTCTATCAA
 TCTCCAGAGAATCCGAGCCCGCCGAGGGCTGCCCCGCCTGGCTGTGTGAGAC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203178 protein sequence
 Red=Cloning site Green=Tags(s)

MAGSREVVAMDCMVGLGPHRESGLARCSLVNVHGAVALYDKFIRPEGEITDYRTRVSGVTPQHVMGATPF
 AVARLEILQLLKGLVVGHDLDKHFQALKEDMSGYTIYDTSTDRLLWREAKLDHCRRVSLRVLSERLLHK
 SIQNSLLGHSSVEDARATMELYQISQIRARRGLPRLAVSD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
 Serotype: AAV-2



[View online »](#)

ACCN:	NM_002201
ORF Size:	543 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002201.4</u> , <u>NP_002192.2</u>
RefSeq Size:	974 bp
RefSeq ORF:	546 bp
Locus ID:	3669
UniProt ID:	<u>Q96AZ6</u>
Cytogenetics:	15q26.1
MW:	20.4 kDa