

## Product datasheet for **RC203192L1V**

### YIPF5 (NM\_030799) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	YIPF5 (NM_030799) Human Tagged ORF Clone Lentiviral Particle
Symbol:	YIPF5
Synonyms:	FinGER5; SB140; SMAP-5; SMAP5; YIP1A
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_030799
ORF Size:	771 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203192).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_030799.6</a>
RefSeq Size:	3212 bp
RefSeq ORF:	774 bp
Locus ID:	81555
UniProt ID:	<a href="#">Q969M3</a>
Cytogenetics:	5q31.3
Domains:	Yip1
Protein Families:	Transmembrane



[View online »](#)

**MW:** 28 kDa

**Gene Summary:** Plays a role in transport between endoplasmic reticulum and Golgi.[UniProtKB/Swiss-Prot Function]