

Product datasheet for **RC203199A1V**

Human COLEC11 (NM_199235) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human COLEC11 (NM_199235) AAV Particle
Tag:	Myc-DDK
Symbol:	COLEC11
Synonyms:	3MC2; CL-11; CL-K1-I; CL-K1-II; CL-K1-IIa; CL-K1-IIb; CLK1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203199 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACGCCTGCTCTGTGCAGATCCTCGTCCCTGGCCTCAAAGGGGATGCGGGAGAGAAGGGAGACAAAGG
CGCCCCCGGACGGCCTGGAAGAGTCGGCCCCACGGGAGAAAAAGCAGTCACAGCCAGTCGTGACAGCTTC
TGACATCTCCAAAAGGAAGTGCACCTCGTCCTTTGTAGAAATGGGATCACAGGGAGACATGGGGGACAAA
GGACAGAAAGGCAGTGTGGGTCGTATGGAAAAATTGGTCCCATTGGCTCTAAAGGTGAGAAAGGAGATT
CCGGTGACATAGGACCCCTGGTCCTAATGGAGAACCAGGCCTCCCATGTGAGTGCAGCCAGCTGCGCAA
GGCCATCGGGGAGATGGACAACCAGGTCTCTCAGCTGACCAGCGAGCTCAAGTTCATCAAGAATGCTGTC
GCCGGTGTGCGCGAGACGGAGAGCAAGATCTACCTGCTGGTGAAGGAGGAGAAGCGCTACGCGGACGCCC
AGCTGTCCTGCCAGGGCCGCGGGGGACGCTGAGCATGCCAAGGACGAGGCTGCCAATGGCCTGATGGC
CGCATACCTGGCGCAAGCCGGCCTGGCCGTGTCTTCATCGGCATCAACGACCTGGAGAAGGAGGGCGCC
TTCGTGTACTCTGACCACTCCCCATGCGGACCTTCAACAAGTGGCGAGCGGTGAGCCCAACAATGCCT
ACGACGAGGAGGACTGCGTGGAGATGGTGGCCTCGGGCGGCTGGAACGACGTGGCCTGCCACACCACCAT
GTACTTCATGTGTGAGTTTGACAAGGAGAACATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203199 protein sequence Red=Cloning site Green=Tags(s) MTPALCRSSSLASKGMRRERRETKAPPDGLEESAPREKKQSQPVVTASDISKRKCTSSFVEMGSQGDMDGK GQKGSVGRHGKIGPIGSKGEKGDSDIGPPGPNGEPGLPCECSQLRKAIGEMDNQVSQLTSELKFIKNAV AGVRETESKIYLLVKEEKRYADAQLSCQGRGGTLSMPKDEAANGLMAAYLAQAGLARVFIGINDLEKEGA FVYSDHSPMRTFNKWRSGEPNNAYDEEDCVEMVASGGWNDVACHTTMYFMCEFDKENM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_199235
ORF Size:	804 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_199235.1</u>
RefSeq Size:	1784 bp
RefSeq ORF:	807 bp
Locus ID:	78989
UniProt ID:	<u>Q9BWP8</u>
Cytogenetics:	2p25.3
MW:	29 kDa