

Product datasheet for **RC203205L3V**

HES6 (NM_018645) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HES6 (NM_018645) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HES6
Synonyms:	bHLHb41; bHLHc23; C-HAIRY1; HES-6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018645
ORF Size:	672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203205).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018645.3
RefSeq Size:	1470 bp
RefSeq ORF:	675 bp
Locus ID:	55502
UniProt ID:	Q96HZ4
Cytogenetics:	2q37.3
Domains:	HLH
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 24.1 kDa

Gene Summary: This gene encodes a member of a subfamily of basic helix-loop-helix transcription repressors that have homology to the Drosophila enhancer of split genes. Members of this gene family regulate cell differentiation in numerous cell types. The protein encoded by this gene functions as a cofactor, interacting with other transcription factors through a tetrapeptide domain in its C-terminus. Alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Dec 2008]