

## Product datasheet for **RC203206L4V**

### ABHD8 (NM\_024527) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ABHD8 (NM_024527) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABHD8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024527
ORF Size:	1317 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203206).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024527.3</a>
RefSeq Size:	2114 bp
RefSeq ORF:	1320 bp
Locus ID:	79575
UniProt ID:	<a href="#">Q96113</a>
Cytogenetics:	19p13.11
Domains:	abhydrolase
MW:	47.3 kDa



[View online »](#)

**Gene Summary:**

This gene is upstream of, and in a head-to-head orientation with the gene for the mitochondrial ribosomal protein L34. The predicted protein contains alpha/beta hydrolase fold and secretory lipase domains. [provided by RefSeq, Jul 2008]