

Product datasheet for **RC203236L2V**

ANP (NPPA) (NM_006172) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ANP (NPPA) (NM_006172) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ANP
Synonyms:	ANF; ANP; ATFB6; ATRST2; CDD; CDD-ANF; CDP; PND
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006172
ORF Size:	459 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203236).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006172.1
RefSeq Size:	858 bp
RefSeq ORF:	456 bp
Locus ID:	4878
UniProt ID:	P01160
Cytogenetics:	1p36.22
Protein Families:	Secreted Protein, Transmembrane
MW:	16.7 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the natriuretic peptide family. Natriuretic peptides are implicated in the control of extracellular fluid volume and electrolyte homeostasis. This protein is synthesized as a large precursor (containing a signal peptide), which is processed to release a peptide from the N-terminus with similarity to vasoactive peptide, cardiodilatin, and another peptide from the C-terminus with natriuretic-diuretic activity. Mutations in this gene have been associated with atrial fibrillation familial type 6. This gene is located adjacent to another member of the natriuretic family of peptides on chromosome 1. [provided by RefSeq, Oct 2015]