

## Product datasheet for **RC203237L2V**

### FMO2 (NM\_001460) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	FMO2 (NM_001460) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FMO2
Synonyms:	FMO1B1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001460
ORF Size:	1413 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203237).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001460.2</a>
RefSeq Size:	5304 bp
RefSeq ORF:	1608 bp
Locus ID:	2327
UniProt ID:	<a href="#">Q99518</a>
Cytogenetics:	1q24.3
Protein Pathways:	Drug metabolism - cytochrome P450
MW:	53.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a flavin-containing monooxygenase family member. It is an NADPH-dependent enzyme that catalyzes the N-oxidation of some primary alkylamines through an N-hydroxylamine intermediate. However, some human populations contain an allele (FMO2\*2A) with a premature stop codon, resulting in a protein that is C-terminally-truncated, has no catalytic activity, and is likely degraded rapidly. This gene is found in a cluster with other related family members on chromosome 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2014]