

## Product datasheet for **RC203245L3V**

### Melanoma Inhibitory Activity (MIA) (NM\_006533) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Melanoma Inhibitory Activity (MIA) (NM_006533) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Melanoma Inhibitory Activity
Synonyms:	CD-RAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006533
ORF Size:	393 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203245).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006533.1</a>
RefSeq Size:	654 bp
RefSeq ORF:	396 bp
Locus ID:	8190
UniProt ID:	<a href="#">Q16674</a>
Cytogenetics:	19q13.2
Domains:	SH3



[View online »](#)

**Protein Families:** Secreted Protein, Transmembrane

**MW:** 14.5 kDa

**Gene Summary:** Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.[UniProtKB/Swiss-Prot Function]