

| | |
|----------------------|---|
| Serotype: | AAV-2 |
| ACCN: | NM_003018 |
| ORF Size: | 591 bp |
| Buffer: | PBS with 0.001% Pluronic F68 |
| Stability: | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| RefSeq: | NM_003018.2 , NP_003009.1 |
| RefSeq Size: | 989 bp |
| RefSeq ORF: | 594 bp |
| Locus ID: | 6440 |
| UniProt ID: | P11686 |
| Cytogenetics: | 8p21.3 |
| MW: | 21.1 kDa |