

Product datasheet for **RC203261L4V**

C15orf15 (RSL24D1) (NM_016304) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C15orf15 (RSL24D1) (NM_016304) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C15orf15
Synonyms:	C15orf15; HRP-L30-iso; L30; RLP24; RPL24; RPL24L; TVAS3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016304
ORF Size:	489 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203261).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016304.2
RefSeq Size:	1487 bp
RefSeq ORF:	492 bp
Locus ID:	51187
UniProt ID:	Q9UHA3
Cytogenetics:	15q21.3
Domains:	Ribosomal_L24e
Protein Pathways:	Ribosome



[View online »](#)

MW: 19.6 kDa

Gene Summary: This gene encodes a protein sharing a low level of sequence similarity with human ribosomal protein L24. Although this gene has been referred to as RPL24, L30, and 60S ribosomal protein L30 isolog in the sequence databases, it is distinct from the human genes officially named RPL24 (which itself has been referred to as ribosomal protein L30) and RPL30. The protein encoded by this gene localizes to the nucleolus and is thought to play a role in the biogenesis of the 60S ribosomal subunit. The precise function of this gene is currently unknown. This gene utilizes alternative polyadenylation signals and has multiple pseudogenes. [provided by RefSeq, Jul 2012]