

## Product datasheet for **RC203271L1V**

### **RAG1AP1 (SLC50A1) (NM\_018845) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RAG1AP1 (SLC50A1) (NM_018845) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAG1AP1
Synonyms:	HsSWEET1; RAG1AP1; SCP; slv; SWEET1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_018845
ORF Size:	663 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203271).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018845.1</a>
RefSeq Size:	1363 bp
RefSeq ORF:	666 bp
Locus ID:	55974
UniProt ID:	<a href="#">Q9BRV3</a>
Cytogenetics:	1q22
Domains:	MtN3_slv
Protein Families:	Transmembrane



[View online »](#)

**MW:** 25 kDa

**Gene Summary:** Mediates sugar transport across membranes. May stimulate V(D)J recombination by the activation of RAG1.[UniProtKB/Swiss-Prot Function]