

Product datasheet for RC203272A1V

Human PEN2 (PSENEN) (NM_172341) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PEN2 (PSENEN) (NM_172341) AAV Particle
Tag: Myc-DDK
Symbol: PEN2
Synonyms: ACNINV2; MDS033; MSTP064; PEN-2; PEN2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203272 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAACCTGGAGCGAGTGTCCAATGAGGAGAAATTGAACCTGTGCCGGAAGTACTACCTGGGGGGTTTG
 CTTTCCTGCCTTTCTCTGGTTGGTCAACATCTTCTGGTTCTTCCGAGAGGCCTTCCTTGCCAGCCTA
 CACAGAACAGAGCCAAATCAAAGGCTATGTCTGGCGCTCAGCTGTGGGCTTCTCTTCTGGGTGATAGTG
 CTCACCTCCTGGATCACCATCTCCAGATCTACCGCCCCGCTGGGGTGCCCTTGGGGAATACCTCTCT
 TCACCATACCCCTGGGCACCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203272 protein sequence
 Red=Cloning site Green=Tags(s)
 MNLERSVNEEKLNLCKRYLGGFAFLPFLWLVNIFWFFREAFVLPAYTEQSIKGYVWRSVGFLEWVIV
 LSWITIFQIYRPRWGALGDYLSFTIPLGTP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_172341
ORF Size: 303 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_172341.1</u>
RefSeq Size:	834 bp
RefSeq ORF:	306 bp
Locus ID:	55851
UniProt ID:	<u>Q9NZ42</u>
Cytogenetics:	19q13.12
MW:	12 kDa