

Product datasheet for **RC203278L1V**

PPT1 (NM_000310) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPT1 (NM_000310) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PPT1
Synonyms:	CLN1; INCL; PPT
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_000310
ORF Size:	918 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203278).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000310.2
RefSeq Size:	2504 bp
RefSeq ORF:	921 bp
Locus ID:	5538
UniProt ID:	P50897
Cytogenetics:	1p34.2
Domains:	Palm_thioest
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Fatty acid elongation in mitochondria, Lysosome, Metabolic pathways

MW: 34.2 kDa

Gene Summary: The protein encoded by this gene is a small glycoprotein involved in the catabolism of lipid-modified proteins during lysosomal degradation. The encoded enzyme removes thioester-linked fatty acyl groups such as palmitate from cysteine residues. Defects in this gene are a cause of infantile neuronal ceroid lipofuscinosis 1 (CLN1, or INCL) and neuronal ceroid lipofuscinosis 4 (CLN4). Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Dec 2008]