

## Product datasheet for **RC203280L1V**

### UTP23 (NM\_032334) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	UTP23 (NM_032334) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UTP23
Synonyms:	C8orf53
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_032334
ORF Size:	747 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203280).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032334.1</a> , <a href="#">NP_115710.1</a>
RefSeq Size:	3693 bp
RefSeq ORF:	750 bp
Locus ID:	84294
UniProt ID:	<a href="#">Q9BRU9</a>
Cytogenetics:	8q24.11
Domains:	DUF652
MW:	28.4 kDa



[View online »](#)

**Gene Summary:** Involved in rRNA-processing and ribosome biogenesis.[UniProtKB/Swiss-Prot Function]