

Product datasheet for RC203286L4

C6orf48 (NM_001040438) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C6orf48 (NM_001040438) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	C6orf48
Synonyms:	D6S57; G8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203286).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001040438
ORF Size:	225 bp



[View online »](#)

