

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for RC203288L4V

## MAGEA3 (NM\_005362) Human Tagged ORF Clone Lentiviral Particle

## **Product data:**

Product Type:	Lentiviral Particles
Product Name:	MAGEA3 (NM_005362) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAGEA3
Synonyms:	CT1.3; HIP8; HYPD; MAGE3; MAGEA6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005362
ORF Size:	942 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203288).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 005362.3</u>
RefSeq Size:	1753 bp
RefSeq ORF:	945 bp
Locus ID:	4102
UniProt ID:	<u>P43357</u>
Cytogenetics:	Xq28
MW:	34.6 kDa



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Gene Summary:This gene is a member of the MAGEA gene family. The members of this family encode<br/>proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the<br/>MAGEA genes show considerable variability, suggesting that the existence of this gene family<br/>enables the same function to be expressed under different transcriptional controls. The<br/>MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in<br/>some hereditary disorders, such as dyskeratosis congenita. [provided by RefSeq, Jul 2008]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US