

Product datasheet for RC203293L2V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

HSD17B2 (NM_002153) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: HSD17B2 (NM_002153) Human Tagged ORF Clone Lentiviral Particle

Symbol: HSD17B2

Synonyms: EDH17B2; HSD17; SDR9C2

Mammalian Cell

Selection:

None

Vector: pLenti-C-mGFP (PS100071)

Tag: mGFP

ACCN: NM_002153 **ORF Size:** 1161 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC203293).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 002153.1</u>

 RefSeq Size:
 1451 bp

 RefSeq ORF:
 1164 bp

 Locus ID:
 3294

 UniProt ID:
 P37059

 Cytogenetics:
 16q23.3

Domains: adh_short

Protein Families: Druggable Genome, Transmembrane



HSD17B2 (NM_002153) Human Tagged ORF Clone Lentiviral Particle - RC203293L2V

Protein Pathways: Androgen and estrogen metabolism, Metabolic pathways

MW: 42.8 kDa

Gene Summary: Capable of catalyzing the interconversion of testosterone and androstenedione, as well as

estradiol and estrone. Also has 20-alpha-HSD activity. Uses NADH while EDH17B3 uses

NADPH.[UniProtKB/Swiss-Prot Function]