

Product datasheet for **RC203295A1V**

Human PBK (NM_018492) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PBK (NM_018492) AAV Particle
Tag:	Myc-DDK
Symbol:	PBK
Synonyms:	CT84; HEL164; Nori-3; SPK; TOPK
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203295 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAAGGGATCAGTAATTTCAAGACACCAAGCAAATTATCAGAAAAAGAAATCTGTATTATGTTCAA
CTCCAATAAATATCCCGCCTCTCCATTTATGCAGAAGCTTGGCTTTGGTACTGGGGTAAATGTGTA
CCTAATGAAAAGATCTCCAAGAGGTTTGTCTCATTCTCCTTGGGCTGTAAAAAGATTAATCCTATATGT
AATGATCATTATCGAAGTGTGTATCAAAAGAGACTAATGGATGAAGCTAAGATTTTGAAGCCCTTCATC
ATCCAAACATTGTTGGTTATCGTGCTTTTACTGAAGCCAATGATGGCAGTCTGTGTCTTGCTATGGAATA
TGGAGGTGAAAAGTCTCTAAATGACTTAATAGAAGAACGATATAAAGCCAGCCAAGATCCTTTTCCAGCA
GCCATAATTTTAAAAGTTGCTTTGAATATGGCAAGAGGGTTAAAGTATCTGCACCAAGAAAAGAACTGC
TTCATGGAGACATAAAGTCTTCAATGTTGTAATTAAGGCGATTTTGAACAATTTAAATCTGTGATGT
AGGAGTCTCTTACCACTGGATGAAAATATGACTGTGACTGACCCTGAGGCTTGTACATTGGCAGAGAG
CCATGGAAACCCAAAGAAGCTGTGGAGGAGAATGGTGTATTACTGACAAGGCAGACATATTTGCCTTTG
GCCTTACTTTGTGGGAAATGATGACTTTATCGATTCCACACATTAATCTTCAAATGATGATGATGATGA
AGATAAACTTTTGTGAAAGTGATTTTGTGATGAAGCATACTATGCAGCGTTGGGAAGTAGGCCACCT
ATTAATATGGAAGAACTGGATGAATCATACCAGAAAGTAATTGAACTCTTCTGTATGCACTAATGAAG
ACCCTAAAGATCGTCCTTCTGCTGCACACATTGTTGAAGCTCTGGAAACAGATGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203295 protein sequence Red=Cloning site Green=Tags(s) MEGISNFKTPSKLSEKKKSVLCSTPTINIPASPFMQKLGFGTGVNVYLMKRSPRGLSHSPWAVKKINPIC NDHYRSVYQKRLMDEAKILKSLHHPNIVGYRAFTEANDGSLCLAMEYGGKSLNDLIEERYKASQDPFPA AIIKVALNMARGLKYLHQEKLLHGDIKSSNVVIKGFETIKICDVGVSLPLDENMTVTDPEACYIGTE PWKPKEAVEENGVIDKADIFAFGLTLWEMMTLSIPHINLSNDDDDKTFDESDFDDEAYYAALGTRPP INMEELDESYQKVIELFSVCTNEDPKDRPSAAHIVEALETDV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_018492
ORF Size:	966 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_018492.2</u>
RefSeq Size:	2127 bp
RefSeq ORF:	969 bp
Locus ID:	55872
UniProt ID:	<u>Q96KB5</u>
Cytogenetics:	8p21.1
MW:	36.1 kDa