

Product datasheet for **RC203295L4V**

PBK (NM_018492) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PBK (NM_018492) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PBK
Synonyms:	CT84; HEL164; Nori-3; SPK; TOPK
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018492
ORF Size:	966 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203295).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018492.2
RefSeq Size:	2127 bp
RefSeq ORF:	969 bp
Locus ID:	55872
UniProt ID:	Q96KB5
Cytogenetics:	8p21.1
Domains:	ptk, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

MW: 36.1 kDa

Gene Summary: This gene encodes a serine/threonine protein kinase related to the dual specific mitogen-activated protein kinase kinase (MAPKK) family. Evidence suggests that mitotic phosphorylation is required for its catalytic activity. The encoded protein may be involved in the activation of lymphoid cells and support testicular functions, with a suggested role in the process of spermatogenesis. Overexpression of this gene has been implicated in tumorigenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]