

Product datasheet for RC203301A1V

Human C1D (NM_173177) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C1D (NM_173177) AAV Particle
Tag: Myc-DDK
Symbol: C1D
Synonyms: hC1D; LRP1; Rrp47; SUN-CoR; SUNCOR
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203301 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGGTGAAGAAATTAATGAAGACTATCCAGTAGAAATTCACGAGTATTTGTCAGCGTTTGAGAATT
 CCATTGGTGCTGTGGATGAGATGCTGAAGACCATGATGTCTGTTTCTAGAAATGAGTTGTTGCAGAAGTT
 GGATCCACTTGAACAAGCAAAAGTGGATTTGGTTTCTGCATACACATTAATCAATGTTTTGGGTTTAT
 TTGGCAACCCAAGGAGTTAATCCTAAGGAACATCCAGTAAACAGGAATTGAAAGAATCAGAGTATATA
 TGAACAGAGTCAAGGAAATAACAGACAAGAAAAAGGCTGGCAAGCTGGACAGAGGTGCAGCTTCAAGATT
 TGTAAAAAATGCCCTCTGGGAACCAAAATCGAAAAATGCATCAAAAGTTGCCAATAAAGGAAAAAGTAA
 AGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203301 protein sequence
 Red=Cloning site Green=Tags(s)

MAGEEINEDYPVEIHEYLSAFENSIGAVDEMLKTMMSVSRNELLQKLDPLEQAKVDLVSAAYTLNSMFWVY
 LATQGVNPKVQELERIRVYMNVRKEITDKKAGKLDRGAASRFVKNALWEPKSKNASKVANKGKSK
 S

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



[View online »](#)

ACCN:	NM_173177
ORF Size:	423 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_173177.1</u>
RefSeq Size:	1194 bp
RefSeq ORF:	426 bp
Locus ID:	10438
UniProt ID:	<u>Q13901</u>
Cytogenetics:	2p14
MW:	16 kDa