

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for RC203307L4V

Cyclophilin A (PPIA) (NM_021130) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cyclophilin A (PPIA) (NM_021130) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cyclophilin A
Synonyms:	CYPA; CYPH; HEL-S-69p
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021130
ORF Size:	495 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203307).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 021130.2</u>
RefSeq Size:	2288 bp
RefSeq ORF:	498 bp
Locus ID:	5478
UniProt ID:	<u>P62937</u>
Cytogenetics:	7p13
Domains:	pro_isomerase
MW:	18 kDa



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Cyclophilin A (PPIA) (NM_021130) Human Tagged ORF Clone Lentiviral Particle – RC203307L4V
Gene Summary:	This gene encodes a member of the peptidyl-prolyl cis-trans isomerase (PPlase) family.
	PPlases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides
	and accelerate the folding of proteins. The encoded protein is a cyclosporin binding-protein
	and may play a role in cyclosporin A-mediated immunosuppression. The protein can also
	interact with several HIV proteins, including p55 gag, Vpr, and capsid protein, and has been

shown to be necessary for the formation of infectious HIV virions. Multiple pseudogenes that

map to different chromosomes have been reported. [provided by RefSeq, Jul 2008]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US