

Product datasheet for **RC203310L4V**

BZW2 (NM_014038) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BZW2 (NM_014038) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BZW2
Synonyms:	HSPC028; MST017; MSTP017
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014038
ORF Size:	1257 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203310).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014038.1
RefSeq Size:	1867 bp
RefSeq ORF:	1260 bp
Locus ID:	28969
UniProt ID:	Q9Y6E2
Cytogenetics:	7p21.1
Domains:	eIF5C
MW:	48.2 kDa



[View online »](#)

Gene Summary: May be involved in neuronal differentiation.[UniProtKB/Swiss-Prot Function]