

Product datasheet for **RC203311L1V**

C10orf7 (CDC123) (NM_006023) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C10orf7 (CDC123) (NM_006023) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C10orf7
Synonyms:	C10orf7; D123
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_006023
ORF Size:	1008 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203311).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006023.1
RefSeq Size:	1566 bp
RefSeq ORF:	1011 bp
Locus ID:	8872
UniProt ID:	O75794
Cytogenetics:	10p14-p13
MW:	39.1 kDa
Gene Summary:	Required for S phase entry of the cell cycle.[UniProtKB/Swiss-Prot Function]



[View online »](#)