

Product datasheet for **RC203314A1V**

Human hnRNP A1 (HNRNPA1) (NM_002136) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human hnRNP A1 (HNRNPA1) (NM_002136) AAV Particle
Tag:	Myc-DDK
Symbol:	hnRNP A1
Synonyms:	ALS19; ALS20; hnRNP-A1; hnRNP A1; HNRPA1; HNRPA1L3; IBMPFD3; UP 1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203314 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCTAAGTCAGAGTCTCCTAAAGAGCCCGAACAGCTGAGGAAGCTCTTCATTGGAGGGTTGAGCTTTG
AAACAACCTGATGAGAGCCTGAGGAGCCATTTTGAGCAATGGGAACGCTCACGGACTGTGTGGTAATGAG
AGATCCAAACACCAAGCGCTCCAGGGGCTTTGGGTTTGTACATATGCCACTGTGGAGGAGGTGGATGCA
GCTATGAATGCAAGGCCACACAAGGTGGATGGAAGAGTTGTGGAACCAAGAGAGCTGTCTCCAGAGAAG
ATTCTCAAAGACCAGGTGCCCACTTAAGTGTGAAAAAGATATTTGTTGGTGGCATTAAAGAAGACACTGA
AGAACATCACCTAAGAGATTATTTGAACAGTATGGAATAATTGAAGTATTGAAATCATGACTGACCGA
GGCAGTGGCAAGAAAAGGGGCTTTGCCTTTGTAACTTTGACGACCATGACTCCGTGGATAAGATTGTCA
TTCAGAAATACCATACTGTGAATGGCCACAACCTGTGAAGTTAGAAAAGCCCTGTCAAAGCAAGAGATGGC
TAGTGCTTCATCCAGCCAAAGAGGTGGAAGTGGTTCTGGAACTTTGGTGGTGGTGGAGGTGGTTTC
GGTGGGAATGACAACCTTCGTCGTGGAGGAACTTCAGTGGTGGTGGTGGCTTTGGTGGCAGCCGTGGT
GTGGTGGATATGGTGGCAGTGGGATGGCTATAATGGATTTGGTAATGATGGAAGCAATTTTGGAGGTGG
TGGAAGCTACAATGATTTTGGGAATTACAACAATCAGTCTTCAAATTTTGGACCATGAAGGGAGGAAAT
TTTGGAGGCAGAAAGCTCTGGCCCTATGGCGGTGGAGGCAATACTTTGCAAAACACGAAACCAAGGTG
GCTATGGCGGTTCCAGCAGCAGCAGTAGCTATGGCAGTGGCAGAAGATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203314 protein sequence Red=Cloning site Green=Tags(s) MSKSESPKEPEQLRKLFIGGLSFETTTDESLRSHFEQWGTLTDCVVMRDPNPKRSRGFGFVTYATVEEVDA AMNARPHKVDGRVVEPKRAVSREDSQRPGAHLTVKKIFVGGIKEDTEHHLRDYFEQYKGIEVIEIMTDR GSGKKRGFAFVTDDHDSVDKIVIQYHTVNGHNCEVRKALSKQEMASASSSQRGRSGSGNFGGGRGGGF GGNDNFGRGGNFSGRGGFGGSRGGGGYGGSGDGYNGFGNDGSNFGGGGSYNDFGNYNQSSNFGPMKGGN FGRSSGPYGGGGQYFAKPRNQGGYGGSSSSSYGSGRRF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002136
ORF Size:	960 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002136.2</u>
RefSeq Size:	1785 bp
RefSeq ORF:	963 bp
Locus ID:	3178
UniProt ID:	<u>P09651</u>
Cytogenetics:	12q13.13
MW:	34.2 kDa