

Product datasheet for **RC203319L3V**

HINT1 (NM_005340) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | HINT1 (NM_005340) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | HINT1 |
| Synonyms: | HINT; NMAN; PKCI-1; PRKCNH1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_005340 |
| ORF Size: | 378 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC203319). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_005340.3 |
| RefSeq Size: | 689 bp |
| RefSeq ORF: | 381 bp |
| Locus ID: | 3094 |
| UniProt ID: | P49773 |
| Cytogenetics: | 5q23.3 |
| Domains: | HIT |
| MW: | 13.8 kDa |


[View online »](#)

Gene Summary:

This gene encodes a protein that hydrolyzes purine nucleotide phosphoramidates substrates, including AMP-morpholidate, AMP-N-alanine methyl ester, AMP-alpha-acetyl lysine methyl ester, and AMP-NH₂. The encoded protein interacts with these substrates via a histidine triad motif. This gene is considered a tumor suppressor gene. In addition, mutations in this gene can cause autosomal recessive neuromyotonia and axonal neuropathy. There are several related pseudogenes on chromosome 7. Several transcript variants have been observed. [provided by RefSeq, Dec 2015]