

Product datasheet for **RC203337A1V**

Human ZFYVE21 (NM_024071) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human ZFYVE21 (NM_024071) AAV Particle
Tag: Myc-DDK
Symbol: ZFYVE21
Synonyms: HCVP7TP1; ZF21
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203337 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCTCCGAGGTGTCCGCGCGCCGCGACGCCAAGAAGCTGGTGCCTCCCGAGCGGCCTGCGCATGG
TGCCCGAACACCGCGCCTTCGGAAGCCCGTTCCGCTGGAGGAGCCGAGTGGGTCCCGACAAGGAGTG
TCGGAGATGTATGCAGTGTGACGCCAAGTTGACTTTCTACCAGAAAGCACCAGTGTGCGCCGTCGCGG
AAGTGCTTCTGCGACAGGTGCTGCAGCCAGAAGGTGCCGCTGCGGCGCATGTGCTTTGTGGACCCCGTGC
GGCAGTGCGCGGAGTGCGCCCTCGTGTCCCTCAAGGAGGCGGAGTTCTACGACAAGCAGCTCAAAGTGCT
CCTGAGCGGAGCCACCTTCTCTGTCACGTTTGAAACTCAGAGAAACCTGAAACTATGACTTGTCTCTT
TCCAATAACCAGAGATACTTGTTTCTGGATGGAGACAGCCACTATGAAATCGAAATTGTACACATTTCCA
CCGTGCAGATCCTCACAGAAGGCTTCCCTCCTGGAGGAGGCAACGCACGGGCCACAGGCATGTTCTGCA
GTATACAGTGCCGGGGACGGAGGGTGTGACCCAGCTGAAGCTGACAGTGGTGGAGGACGTGACTGTGGGC
AGGAGGCAGGCGGTGGCGTGGCTAGTGGCCATGCACAAGGCTGCCAAGCTCCTCTATGAATCTCGGGACC
AG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203337 protein sequence Red=Cloning site Green=Tags(s) MSSEVSARRDAKKLVSPSGLRMVPEHRAFGSPFGLLEPQWVPDKECRRCMQCDAKFDLTRKHHCRRCG KCFCDRCCSQKVPLRRMCFVDPVRQCAECALVSLKEAEFYDKQLKVLLSGATFLVTFGNSEKPEMTTCRL SNNQRYLFLDGDShYEIEIVHISTVQILTEGFPPGGGNARATGMFLQYTPGTEGVTQLKLTVEDVTVG RRQAVAWLVAMHKAAKLLYESRDQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024071
ORF Size:	702 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024071.2</u>
RefSeq Size:	1470 bp
RefSeq ORF:	705 bp
Locus ID:	79038
UniProt ID:	<u>Q9BQ24</u>
Cytogenetics:	14q32.33
MW:	26.5 kDa