

Product datasheet for **RC203343L3V**

RPL3 (NM_000967) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | RPL3 (NM_000967) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | RPL3 |
| Synonyms: | ASC-1; L3; TARBP-B |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_000967 |
| ORF Size: | 1209 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC203343). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_000967.3 |
| RefSeq Size: | 1348 bp |
| RefSeq ORF: | 1212 bp |
| Locus ID: | 6122 |
| UniProt ID: | P39023 |
| Cytogenetics: | 22q13.1 |
| Domains: | Ribosomal_L3 |
| Protein Pathways: | Ribosome |



[View online »](#)

MW: 46.1 kDa

Gene Summary: Ribosomes, the complexes that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L3P family of ribosomal proteins and it is located in the cytoplasm. The protein can bind to the HIV-1 TAR mRNA, and it has been suggested that the protein contributes to tat-mediated transactivation. This gene is co-transcribed with several small nucleolar RNA genes, which are located in several of this gene's introns. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]