

Product datasheet for **RC203347A1V**

Human MNF1 (UQCC2) (NM_032340) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human MNF1 (UQCC2) (NM_032340) AAV Particle

Tag: Myc-DDK

Symbol: MNF1

Synonyms: bA6B20.2; C6orf125; C6orf126; Cbp6; M19; MC3DN7; MNF1

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC203347 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGCCAGCCGGTACCGCGTTTTCTTAAGCTCTGTGAGGAATGGCCAGTGGACGAGACCAACGGG
GCCGGGACTTGGGCGTTACCTGCGACAGCGGTAGCACAGGCCTTTCGGGAGGGAGAGAATACCCAGGT
TGCAGAGCCTGAGGCCTGTGATCAGATGTACGAGAGCTTAGCGCGACTCCATTCAAACACTACAACAC
AAGTACCCTCGCCCCAGAGACACCAGCTTCAGTGGCCTGTCGTTGGAAGAGTACAAGCTGATCCTGTCCA
CAGACACCTTGAAGAGCTTAAGGAAATAGATAAAGGCATGTGGAAGAACTGCAGGAGAAGTTTGCCCC
CAAGGGTCCTGAGGAGGATCATAAGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203347 protein sequence
Red=Cloning site Green=Tags(s)

MAASRYRRFLKLCEEWPVDETKRGRDLGAYLRQVAQAFREAGENTQVAEPEACDQMYESLARLHSNYYKH
KYPRPRDTSFSGLSLEEKYKILSTDTLEELKEIDKGMWKKLQEKFAFKGPEEDHKA

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_032340



[View online »](#)

ORF Size:	378 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032340.1</u>
RefSeq Size:	1355 bp
RefSeq ORF:	381 bp
Locus ID:	84300
UniProt ID:	<u>Q9BRT2</u>
Cytogenetics:	6p21.31
MW:	14.9 kDa