

## Product datasheet for **RC203367L4V**

### **C2orf43 (LDAH) (NM\_021925) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	C2orf43 (LDAH) (NM_021925) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C2orf43
Synonyms:	C2orf43; hLDAH
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021925
ORF Size:	975 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203367).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021925.1</a>
RefSeq Size:	3963 bp
RefSeq ORF:	978 bp
Locus ID:	60526
UniProt ID:	<a href="#">Q9H6V9</a>
Cytogenetics:	2p24.1
Protein Families:	Druggable Genome
MW:	37.3 kDa



[View online »](#)

**Gene Summary:**

Serine lipid hydrolase associated with lipid droplets. Highly expressed in macrophage-rich areas in atherosclerotic lesions, suggesting that it could promote cholesterol ester turnover in macrophages.[UniProtKB/Swiss-Prot Function]