

Product datasheet for **RC203380L4V**

DUS3L (NM_020175) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DUS3L (NM_020175) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DUS3L
Synonyms:	DUS3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020175
ORF Size:	1950 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203380).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020175.1 , NP_064560.1
RefSeq Size:	2127 bp
RefSeq ORF:	1953 bp
Locus ID:	56931
UniProt ID:	Q96G46
Cytogenetics:	19p13.3
Domains:	zf-CCCH, Dus
MW:	72.5 kDa



[View online »](#)

Gene Summary:

Catalyzes the synthesis of dihydrouridine, a modified base found in the D-loop of most tRNAs.[UniProtKB/Swiss-Prot Function]