

Product datasheet for **RC203392A1V**

Human CHODL (NM_024944) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CHODL (NM_024944) AAV Particle
Tag:	Myc-DDK
Symbol:	CHODL
Synonyms:	C21orf68; MT75; PRED12
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203392 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGCCGCGTGGTCTCGCTGCTGCTGGGCGCCGCGCTGCTCTGCGGCCACGGAGCCTTCTGCCGCCGCG
TGGTCAGCGGCCAAAAGGTGTGTTTTGCTGACTTCAAGCATCCCTGCTACAAATGGCCTACTTCCATGA
ACTGTCCAGCCGAGTGAGCTTTCAGGAGGCACGCTGGCTTGAGAGTGAGGGAGGAGTCTCTCAGC
CTTGAGAATGAAGCAGAACAGAAGTTAATAGAGAGCATGTTGCAAAACCTGACAAACCCGGGACAGGGA
TTTCTGATGGTGATTTCTGGATAGGGCTTTGGAGGAATGGAGATGGGCAACATCTGGTGCCTGCCAGA
TCTCTACCACTGGTCTGATGGAAGCAATTCACGATACCGAAACTGGTACACAGATGAACCTTCTGCGGA
AGTGAAAAGTGTGTTGTGATGTATACCAACCAACTGCCAATCCTGGCCTTGGGGTCCCTACCTTTACC
AGTGGAATGATGACAGGTGTAACATGAAGCACAATTATTTGCAAGTATGAACCAGAGATTAATCCAAC
AGCCCCTGTAGAAAAGCCTTATCTTACAAATCAACCAGGAGACACCCATCAGAATGTGGTTGTTACTGAA
GCAGGTATAATTCCCAATCTAATTTATGTTGTTATACCAACAATACCCCTGCTCTTACTGATACTGGTTG
CTTTTGAACCTGTTGTTTCCAGATGCTGCATAAAAGTAAAGGAAGAACAAAAGTATGCCAAACCATGTC
TACACTGTGGATTTCAAAGAGTACCAGAAAAGAAAGTGGCATGGAAGTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203392 protein sequence Red=Cloning site Green=Tags(s) MSRVVSLLLGAALLCGHGAFRRVVSQGKVCFAADFKHPCYKMAFHELSSRVSFQEARLACESEGGVLLS LENEAEQKLIESMLQNLTKPGTGISDGDWIGLWRNGDGTSGACPDLYQWSDGSNSQYRNWYTDEPSCG SEKCVVMYHQPTANPGLGGPYLYQWDDRCNMKHNYICKYEPEINPTAPVEKPYLTNQPGDTHQNVVTE AGIIPNLIYVVIPTIPLLLLILVAFGTCCFQMLHKSCKGRKTSPNQSTLWISKSTRKESGMEV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024944
ORF Size:	819 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024944.2</u>
RefSeq Size:	2625 bp
RefSeq ORF:	822 bp
Locus ID:	140578
UniProt ID:	<u>Q9H9P2</u>
Cytogenetics:	21q21.1
MW:	30.4 kDa