

Product datasheet for **RC203407L1V**

Artemis (DCLRE1C) (NM_022487) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Artemis (DCLRE1C) (NM_022487) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Artemis
Synonyms:	A-SCID; DCLREC1C; RS-SCID; SCIDA; SNM1C
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_022487
ORF Size:	1731 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203407).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022487.2
RefSeq Size:	6242 bp
RefSeq ORF:	1734 bp
Locus ID:	64421
UniProt ID:	Q96SD1
Cytogenetics:	10p13
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Non-homologous end-joining, Primary immunodeficiency



[View online »](#)

MW: 65.3 kDa

Gene Summary: This gene encodes a nuclear protein that is involved in V(D)J recombination and DNA repair. The encoded protein has single-strand-specific 5'-3' exonuclease activity; it also exhibits endonuclease activity on 5' and 3' overhangs and hairpins. The protein also functions in the regulation of the cell cycle in response to DNA damage. Mutations in this gene can cause Athabaskan-type severe combined immunodeficiency (SCIDA) and Omenn syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]