

Product datasheet for RC203419A1V

Human Destrin (DSTN) (NM_006870) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Destrin (DSTN) (NM_006870) AAV Particle
Tag: Myc-DDK
Symbol: Destrin
Synonyms: ACTDP; ADF; bA462D18.2; HEL32
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203419 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGCCTCAGGAGTGCAAGTAGCTGATGAAGTATGTCGCATTTTTATGACATGAAAGTTCGTAATGCT
 CCACACCAGAAGAAATCAAGAAAAGAAAGAGGCTGTCATTTTTGTCTCAGTGCAGACAAAAGTGCAT
 CATTGTAGAAGAAGGCAAGAGATCTTGGTTGGAGATGTTGGTGAACCATAACTGATCCTTTCAAGCAT
 TTTGTGGGAATGCTTCTGAAAAAGATTGTCGCTATGCTTTGTATGATGCAAGCTTTGAAACAAAAGAT
 CCAGAAAAGAAGAGTTGATGTTTTTTTTGTGGCACCAGAACTAGCACCTCTGAAAAGTAAATGATCTA
 TGCAAGCTCCAAGGATGCAATTAAGAAAGAAATTTCAAGGCATAAAACATGAATGTCAAGCAAAATGGACCA
 GAAGATCTCAATCGGGCTTGATTGCTGAAAAGTTAGGTGGATCCTTAATTGTAGCCTTTGAAGGATGCC
 CTGTG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203419 protein sequence
 Red=Cloning site Green=Tags(s)

MASGVQVADEVCRIFYDMKVRKSTPEEIKRKKAVIFCLSDKKCIIVEEGKEILVGDVGVITDPFKH
 FVGMLEPKDCRYALYDASFETKESRKEELMFFLWAPELAPLKSMMIYASSKDAIKKKFQGIKHECQANGP
 EDLNRACIAEKLGGSLIVAFEGCPV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_006870
ORF Size:	495 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006870.3
RefSeq Size:	1614 bp
RefSeq ORF:	498 bp
Locus ID:	11034
UniProt ID:	P60981
Cytogenetics:	20p12.1
MW:	18.5 kDa