

Product datasheet for **RC203419L3V**

Destrin (DSTN) (NM_006870) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Destrin (DSTN) (NM_006870) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Destrin
Synonyms:	ACTDP; ADF; bA462D18.2; HEL32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006870
ORF Size:	495 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203419).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006870.3
RefSeq Size:	1614 bp
RefSeq ORF:	498 bp
Locus ID:	11034
UniProt ID:	P60981
Cytogenetics:	20p12.1
Domains:	ADF
MW:	18.5 kDa



[View online »](#)

Gene Summary:

The product of this gene belongs to the actin-binding proteins ADF family. This family of proteins is responsible for enhancing the turnover rate of actin in vivo. This gene encodes the actin depolymerizing protein that severs actin filaments (F-actin) and binds to actin monomers (G-actin). Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]