

Product datasheet for **RC203427L4V**

DCUN1D1 (NM_020640) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DCUN1D1 (NM_020640) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DCUN1D1
Synonyms:	DCNL1; DCUN1L1; RP42; SCCRO; SCRO; Tes3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020640
ORF Size:	777 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203427).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020640.2
RefSeq Size:	3192 bp
RefSeq ORF:	780 bp
Locus ID:	54165
UniProt ID:	Q96GG9
Cytogenetics:	3q26.33
Domains:	DUF298
Protein Families:	Druggable Genome



[View online »](#)

MW: 30.1 kDa

Gene Summary: Part of an E3 ubiquitin ligase complex for neddylation. Promotes neddylation of cullin components of E3 cullin-RING ubiquitin ligase complexes. Acts by binding to cullin-RBX1 complexes in the cytoplasm and promoting their nuclear translocation, enhancing recruitment of E2-NEDD8 (UBE2M-NEDD8) thioester to the complex, and optimizing the orientation of proteins in the complex to allow efficient transfer of NEDD8 from the E2 to the cullin substrates. Involved in the release of inhibitory effects of CAND1 on cullin-RING ligase E3 complex assembly and activity (PubMed:25349211, PubMed:28581483). Acts also as an oncogene facilitating malignant transformation and carcinogenic progression (By similarity). [UniProtKB/Swiss-Prot Function]