

Product datasheet for **RC203429L4V**

Sortilin (SORT1) (NM_002959) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sortilin (SORT1) (NM_002959) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SORT1
Synonyms:	Gp95; LDLCQ6; NT3; NTR3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002959
ORF Size:	2493 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203429).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002959.4 , NP_002950.3
RefSeq Size:	7018 bp
RefSeq ORF:	2496 bp
Locus ID:	6272
UniProt ID:	Q99523
Cytogenetics:	1p21.3-p13.1
Domains:	BNR, VPS10
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Protein Pathways: Lysosome, Neurotrophin signaling pathway

MW: 92.07 kDa

Gene Summary: This gene encodes a member of the VPS10-related sortilin family of proteins. The encoded preproprotein is proteolytically processed by furin to generate the mature receptor. This receptor plays a role in the trafficking of different proteins to either the cell surface, or subcellular compartments such as lysosomes and endosomes. Expression levels of this gene may influence the risk of myocardial infarction in human patients. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2015]